


# Disposable Surgical Masks





<b>BMR</b>	<b>4</b>	<b>Klear Vision™</b>
<b>FLUID PROTECTION</b>		
<b>Maximum Fluid Protection Masks with Eyeshield</b> <ul style="list-style-type: none"> <li>Non-fogging wraparound optical quality eyeshield assures maximum visibility and facial protection from impacting fluids.</li> <li>A two-ply filter in each facemask provides maximum resistance to fluid impact.</li> <li>Cellulose innerfacing for skin comfort.</li> <li>Comfortable anti-fog foam strip prevents fogging of eyewear.</li> <li>Aluminum nose wire for secure fit.</li> </ul>		
Item #	Description	Pkg.
3505PA	Max Fluid Protection Anti-fog Tie & Eyeshield	25/Box, 100/cs
3508PA	Max Fluid Protection Anti-fog Ear Loop & Eyeshield	25/Box, 100/cs
<b>Maximum Fluid Protection Masks</b> <ul style="list-style-type: none"> <li>A two-ply filter in each facemask provides maximum resistance to fluid impact.</li> <li>Cellulose innerfacing for skin comfort.</li> <li>Comfortable anti-fog foam strip prevents fogging of eyewear*.</li> <li>Aluminum nose wire for secure fit.</li> </ul>		
Item #	Description	Pkg.
3503P	Max Fluid Protection Tie	50/Box, 300/cs
3503PA	Max Fluid Protection Anti-Fog Tie*	50/Box, 300/cs
3501PA	Max Fluid Protection Anti-Fog Ear Loop*	50/Box, 300/cs
<b>Fluid Protection Masks with Eyeshield</b> <ul style="list-style-type: none"> <li>Non-fogging wraparound optical quality eyeshield assures maximum visibility and facial protection from impacting fluids.</li> <li>An extra protection filter provides excellent resistance to fluid impact.</li> <li>Cellulose innerfacing for skin comfort.</li> <li>Comfortable anti-fog foam strip prevents fogging of eyewear.</li> <li>Aluminum nose wire for secure fit.</li> </ul>		
Item #	Description	Pkg.
2505BA	Fluid Protection Anti-fog Tie & Eyeshield	25/Box, 100/cs
2508BA	Fluid Protection Anti-fog Ear Loop & Eyeshield	25/Box, 100/cs
<b>Fluid Protection Masks</b> <ul style="list-style-type: none"> <li>An extra protection filter provides excellent resistance to fluid impact.</li> <li>Cellulose innerfacing for skin comfort.</li> <li>Comfortable anti-fog foam strip prevents fogging of eyewear.</li> <li>Aluminum nose wire for secure fit.</li> </ul>		
Item #	Description	Pkg.
2503BA	Fluid Protection Anti-fog Tie	50/Box, 300/cs
2501BA	Fluid Protection Anti-Fog Ear Loop	50/Box, 300/cs




<b>BMR</b>	<b>5</b>	<b>Klear Vision™</b>						
<p><b>Anti-Fog</b> You have choices</p>								
<p><b>Anti-fog Adhesive Tie Mask</b></p> <ul style="list-style-type: none"> <li>· Medical grade tape strip to prevent fogging of eyewear.</li> <li>· Cellulose innerfacing for skin comfort.</li> <li>· Aluminum nose wire for secure fit.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1624G</td> <td>Anti-fog Adhesive Tie</td> <td>50/Box,300/cs</td> </tr> </tbody> </table>			Item #	Description	Pkg.	1624G	Anti-fog Adhesive Tie	50/Box,300/cs
Item #	Description	Pkg.						
1624G	Anti-fog Adhesive Tie	50/Box,300/cs						
<p><b>Anti-fog Film Tie Mask</b></p> <ul style="list-style-type: none"> <li>· EMA film strip to prevent fogging of eyewear.</li> <li>· Cellulose innerfacing for skin comfort.</li> <li>· Aluminum nose wire for secure fit.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1504G</td> <td>Anti-fog Film Tie</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1504G	Anti-fog Film Tie	50/Box, 300/cs	
Item #	Description	Pkg.						
1504G	Anti-fog Film Tie	50/Box, 300/cs						
<p><b>Anti-fog Foam Tie Mask</b></p> <ul style="list-style-type: none"> <li>· Soft foam strip to prevent fogging of eyewear.</li> <li>· Cellulose innerfacing for skin comfort.</li> <li>· Aluminum nose wire for secure fit.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1514G</td> <td>Anti-fog Foam Tie</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1514G	Anti-fog Foam Tie	50/Box, 300/cs	
Item #	Description	Pkg.						
1514G	Anti-fog Foam Tie	50/Box, 300/cs						

<b>BMR</b>	<b>6</b>	<b>Klear Vision™</b>							
<b>Specialty</b>									
<p><b>Chamber Mask</b></p> <ul style="list-style-type: none"> <li>· Largest breathing chamber.</li> <li>· Cellulose innerfacing for skin comfort.</li> <li>· Full width nosepiece for proper fit.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1606G</td> <td>Chamber Style Tie</td> <td>50/Box, 300</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1606G	Chamber Style Tie	50/Box, 300		
Item #	Description	Pkg.							
1606G	Chamber Style Tie	50/Box, 300							
<p><b>Sensitive Skin</b></p> <ul style="list-style-type: none"> <li>· Contains no inks or colorants that could irritate the skin.</li> <li>· Cellulose innerfacing for skin comfort.</li> <li>· Aluminum nose wire for secure fit.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1502W</td> <td>Sensitive Skin Tie</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1502W	Sensitive Skin Tie	50/Box, 300/cs		
Item #	Description	Pkg.							
1502W	Sensitive Skin Tie	50/Box, 300/cs							
<p><b>Laser Mask</b></p> <ul style="list-style-type: none"> <li>· The extra protection filter is rated at 0.1 micron particle size.</li> <li>· Cellulose innerfacing for skin comfort.</li> <li>· Aluminum nose wire for secure fit.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>2510L</td> <td>Laser Tie</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	2510L	Laser Tie	50/Box, 300/cs		
Item #	Description	Pkg.							
2510L	Laser Tie	50/Box, 300/cs							

<b>BMR</b>	7	<b>Klear Vision™</b>									
<b>Surgical</b> The standard for the O.R.											
<b>Standard Surgical Tie Mask</b> <ul style="list-style-type: none"> <li>Cellulose innerfacing for skin comfort.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1503B</td> <td>Standard Surgical Tie</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1503B	Standard Surgical Tie	50/Box, 300/cs				
Item #	Description	Pkg.									
1503B	Standard Surgical Tie	50/Box, 300/cs									
<b>Lite Surgical Tie Mask</b> <ul style="list-style-type: none"> <li>Our Lightest.</li> <li>Cellulose innerfacing for skin comfort.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1403B</td> <td>Lite Surgical Tie Blue</td> <td>50/Box, 300/cs</td> </tr> <tr> <td>1403G</td> <td>Lite Surgical Tie Green</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1403B	Lite Surgical Tie Blue	50/Box, 300/cs	1403G	Lite Surgical Tie Green	50/Box, 300/cs	
Item #	Description	Pkg.									
1403B	Lite Surgical Tie Blue	50/Box, 300/cs									
1403G	Lite Surgical Tie Green	50/Box, 300/cs									
<b>Soft Surgical Tie Mask</b> <ul style="list-style-type: none"> <li>Our Softest.</li> <li>Bi-Component innerfacing for maximum softness.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1303B</td> <td>Soft Surgical Tie</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1303B	Soft Surgical Tie	50/Box, 300/cs				
Item #	Description	Pkg.									
1303B	Soft Surgical Tie	50/Box, 300/cs									
<b>Cone Style Mask</b> <ul style="list-style-type: none"> <li>Easy to Don.</li> <li>Non-collapsing Breathing Chamber.</li> <li>Soft and Comfortable.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1550</td> <td>Cone Style Mask</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1550	Cone Style Mask	50/Box, 300/cs				
Item #	Description	Pkg.									
1550	Cone Style Mask	50/Box, 300/cs									

<b>BMR</b>	<b>8</b>	<b>Klear Vision™</b>									
<b>Procedure Masks</b> The Standard Ear Loop											
<b>Procedure Mask</b> <ul style="list-style-type: none"> <li>Tissue innerfacing for comfort.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1001B</td> <td>Procedure Mask</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>			Item #	Description	Pkg.	1001B	Procedure Mask	50/Box, 300/cs			
Item #	Description	Pkg.									
1001B	Procedure Mask	50/Box, 300/cs									
<b>Lite Procedure Mask</b> <ul style="list-style-type: none"> <li>Our Lightest.</li> <li>Cellulose innerfacing for skin comfort.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1401B</td> <td>Procedure Mask</td> <td>50/Box, 300/cs</td> </tr> <tr> <td>1401G</td> <td>Procedure Mask Green</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1401B	Procedure Mask	50/Box, 300/cs	1401G	Procedure Mask Green	50/Box, 300/cs	
Item #	Description	Pkg.									
1401B	Procedure Mask	50/Box, 300/cs									
1401G	Procedure Mask Green	50/Box, 300/cs									
<b>Procedure Mask—Mint Scented</b> <ul style="list-style-type: none"> <li>Cellulose innerfacing for skin comfort.</li> <li>Mint Scented for comfort.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1401GM</td> <td>Procedure Mask</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1401GM	Procedure Mask	50/Box, 300/cs				
Item #	Description	Pkg.									
1401GM	Procedure Mask	50/Box, 300/cs									
<b>Isolation Mask</b> <ul style="list-style-type: none"> <li>Tissue innerfacing for comfort.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>1001Y</td> <td>Procedure Mask</td> <td>50/Box, 300/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	1001Y	Procedure Mask	50/Box, 300/cs				
Item #	Description	Pkg.									
1001Y	Procedure Mask	50/Box, 300/cs									

Description	Kimberly-Clark	Berkley Medical
Max Fluid Protection Anti-fog Tie & Eye Shield	48247	3505PA
Max Fluid Protection Anti-fog Ear Loop & Eye Shield	47137	3508PA
Max Fluid Protection Tie	48208	3503P
Max Fluid Protection Anti-Fog Tie	48207	3503PA
Max Fluid Protection Anti-Fog Ear Loop	47107	3501PA
Fluid Protection Anti-fog Tie & Eye Shield	62114	2505BA
Fluid Protection Anti-fog Ear Loop & Eye Shield	62116	2508BA
Fluid Protection Anti-fog Tie	62113	2503BA
Fluid Protection Anti-fog Ear Loop	62115	2501BA
Standard Ear Loop		1501B
Sensitive Skin Tie	49310	1502W
Standard Surgical Tie	48201	1503B
Anti-fog Film Tie	49230	1504G
Anti-fog Foam Tie	49214	1514G
Anti-fog Adhesive Tie	49215	1624G
Laser Tie	47625	2510L
Lite Surgical Tie	48100	1403B
Lite Surgical Tie	48105	1403G
Soft Surgical Tie	48390	1303B
Cone Style	47200	1550
Lite Procedure Ear Loop	62356	1401B
Lite Procedure Ear Loop	62357	1401G
Procedure Ear Loop—Mint Scented		1401GM
Procedure Ear Loop	47080	1001B
Isolation Ear Loop	47117	1001Y
N95 Cone Respirator	46727	6001
N95 Flat Fold Respirator		6101

<b>BMR</b>	<b>9</b>	<b>Klear Vision™</b>						
N95								
<p>N95</p> <ul style="list-style-type: none"> <li>· Comfortable innerfacing</li> <li>· Contoured for facial seal.</li> <li>· Aluminum nosepiece for proper fit.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>6001</td> <td>N95 Cone Respirator</td> <td>20/Box, 120/cs</td> </tr> </tbody> </table>		Item #	Description	Pkg.	6001	N95 Cone Respirator	20/Box, 120/cs	
Item #	Description	Pkg.						
6001	N95 Cone Respirator	20/Box, 120/cs						
<p>N95</p> <ul style="list-style-type: none"> <li>· Comfortable innerfacing</li> <li>· Contoured for facial seal.</li> <li>· Aluminum nosepiece for proper fit.</li> </ul> <table border="1"> <thead> <tr> <th>Item #</th> <th>Description</th> <th>Pkg.</th> </tr> </thead> <tbody> <tr> <td>6101</td> <td>N95 Flat Fold Respirator</td> <td>20/Box, 120/cs</td> </tr> </tbody> </table>			Item #	Description	Pkg.	6101	N95 Flat Fold Respirator	20/Box, 120/cs
Item #	Description	Pkg.						
6101	N95 Flat Fold Respirator	20/Box, 120/cs						



## GLOSSARY

### **BFE (BACTERIAL FILTRATION EFFICIENCY %)**

The efficiency at which the facemask barrier prevents aerosolized bacteria penetration at a given flow rate.

### **PFE (PARTICLE FILTRATION EFFICIENCY %)**

The efficiency at which the facemask barrier prevents aerosolized particle penetration at a given flow rate.

### **DELTA P ( $\Delta P$ )**

The pressure drop across the facemask barrier that occurs when simulating human breathing and is measured in  $\text{mm H}_2\text{O} / \text{CM}^2$  at a given flow rate

### **FLUID PROTECTION**

The resistance of the facemask barrier to fluid impact penetration as rated at pressures and flow rates simulating the human body. The ASTM F 1862 test has 3 levels of protection; 80mm, 120mm AND 160mm.

### **PROCEDURE FACEMASK**

A facemask with ear loops.

### **SURGICAL FACEMASK**

A facemask with tie strings.

### **N95 FACEMASK**

A facemask that provides a leak-proof facial fit and meets the barrier requirements of NIOSH and FDA.

### **ANTI-FOG ADHESIVE TAPE**

A heat-activated medical grade adhesive tape that forms a vapor seal that prevents warm, humid exhaled breath from fogging eye wear.

### **ANTI-FOG ADHESIVE FOAM**

A medical grade foam that conforms to the wearer to minimize warm, humid exhaled breath from fogging eye wear.

### **ANTI-FOG FILM**

An EMA film that forms a vapor seal that prevents warm, humid exhaled breath from fogging eye wear.

### **EYESHIELD**

An optical quality plastic visor that is attached to fluid protection facemasks to prevent fluids from contacting the wearer's eyes.

### **ANTI-GLARE STRIP**

A dark colored strip at the top of the facemask that removes glare from the eyeshield.

### **BARRIER**

The material used to cover the face. The barrier most commonly consists of 3 plies: innerfacing, filter and outerfacing. The innerfacing must comply with ISO 10993-1 and be comfortable to the skin. The filter is generally polypropylene melt-blown and provides the greatest part of the protection. The outerfacing provides additional protection, covers the filter to prevent dusting and provides color and style. The barrier must meet FDA requirements of 95% to 100% B.F.E. at 3.0 microns and be less than 5.0 Delta P for breathing. Barriers that claim fluid protection compliance with ASTM 2100 must be tested to ASTM F.1862 .